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Beispiele

PCR-Step	Glass	Polycarb.	Glass	Polycarb.	Glass	Polycarb.
Denaturierung	0 sec.	7 sec.	5 sec	7-10 sec	>7 sec	>7 sec
Annealing	0 sec.	7 sec	5 sec	7-10 sec	>7 sec	>7 sec

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Genaxxon BioScience

Polycarbonate Capillaries User Guide

Product	Cat#	Package size
Polycarbonate Capillaries for use in capillary PCR	I2250.0960	960 capillaries
Starter Kit with polycarbonate capillaries	I2250.0961	672 capillaries + 1 carousel
PEEK-Carousel for use in LightCycler machines	I2250.0001	1 carousel

Sample preparation / PCR Mixture preparation

Samples for PCR can be handled and prepared according to the existing protocol at customer site. No changes have to be made using polycarbonate capillaries instead the original glass capillaries from Hofmann La-Roche.

No changes have to be made preparing the PCR Mix.

It is possible to use our polycarbonate capillaries for PCR volumes from 10 µl to 50 µl without changes in the PCR programme (for details see point PCR programme).

Handling capillaries and caps

For better/easier handling of the capillaries and caps we provide a "Transfer-Pin" for taking the capillaries and the caps out of the rack-box. The Transfer-Pin is part of our "Starter-Kit". If you are used to work without a "Transfer-Pin", you can also work without it, but use protective gloves to protect you PCR-reaction from contamination, by for example DNAses.

For sealing the capillaries with the caps, please use the transfer pin, to transfer caps from the rack-box to the polycarbonate capillary. Ensure that each capillary is closed tightly by checking it visually (the smaller part of the caps must be completely inserted into the plastic capillary). Normally you can hear 2 times a "click", if the cap fits well.

The LightCycler Carousel(s)

Genaxxon offers 2 versions of LightCycler carousel that fit to LightCycler Versions 1.2 and 1.5 or to Version 2.0 respectively. Same with the centrifuge rotor for the LightCycler 1.2 / 1.5 or 2.0 carousel. It is not possible to use the carousels in the wrong Version of a LightCycler or to use the wrong carousel in the wrong centrifuge rotor.

Please inform Genaxxon about the LightCycler version available in your laboratory.

Performing PCR / PCR programme

Our polycarbonate capillaries have a bigger diameter, compared to the 20 µl glass capillaries from Roche.

To guarantee effective heat transfer from the outside of the capillaries to the middle of the PCR reaction, despite the bigger diameter, it might be necessary to increase the "hold-times" at Denaturing and Annealing.

For this reason we recommend strongly to program your LightCycler® with at least 7-10 seconds at 95°C for the denaturing step and 7 seconds for the annealing step. Too short denaturing and annealing times will lead to a complete failure of your PCR, or at least to much higher CT-times.

If you are already using a programme with both times higher than 7 seconds, no changes have to be made normally. Nevertheless it might be useful to increase the existing denaturing and/or annealing times by 3 to 5 seconds each.

As the easiness of denaturing DNA depends on the GC-content of your DNA, it might be better to prolong denaturing up to 15 seconds (only in case of **NO** PCR signal with a denaturing time of 10 seconds).

Example

PCR-Step	Glass	Polycarb.	Glass	Polycarb.	Glass	Polycarb.
Denaturing	0 sec.	7 sec.	5 sec	7-10 sec	>7 sec	>7 sec
Annealing	0 sec.	7 sec	5 sec	7-10 sec	>7 sec	>7 sec

Sample recovery after PCR

Some customer want to analyse the PCR reaction on a gel. For that purpose the cap can be taken from the capillary with the Transfer-Pin. The recovery of the sample is principally possible using pipettes with a 10 µl tip. In some cases the pipet tips are too broad (Sarstedt, AHN and other brands). In this case you can cut the capillaries with scissors and spin the solution into a clean reaction tube.

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