

# IX. Consumables

## Product Overview

### Gloves

- Latex
- Nitril
- Vinyl
- Miscellaneous

### Protective Masks

### Tubes

- Cryo Tubes
- Falcon Type Tubes
- Reaction Tubes
- PCR Tubes

### Pipet Tips & Serological Pipets

- Pipet Tips
- Filter Tips
- Dispenser Tips
- Serological Pipets

### Racks & Boxes

- Racks
- Freezer Boxes
- Card-board Freezer Boxes

### Tissue Culture

## Product Overview

### Gloves

#### Latex

Prod.No.	Product	Pack-size	Price in €
M6138.10xx	Blossom Latex Aloe Vera Plus	10 x 100 Stück	80.27
M6139.10xx	Latex Superior with mint	10 x 100 Stück	49.70
M6142.005S	Prima Duo	25 Paare	10.39
M6142.005M	Prima Duo	25 Paare	10.39
M6142.005L	Prima Duo	25 Paare	10.39
M6143.10xx	Prima Plus Blau	10 x 100 Stück	51.90
M6144.10xx	Prima Plus Natur	10 x 100 Stück	51.90
M6376.04xx	Sterile Operationshandschuhe (rechts/links)	4 x 50 Paar	129.20
M6376.20xx	Sterile Operationshandschuhe (rechts/links)	20 x 50 Paar	543.77
M6374.06xx	Sterile Operationshandschuhe (rechts/links)	6 x 50 Paar	191.11
M6374.30xx	Sterile Operationshandschuhe (rechts/links)	30 x 50 Paar	813.98
M6375.06xx	Sterile Operationshandschuhe (rechts/links)	6 x 50 Paar	183.36
M6375.30xx	Sterile Operationshandschuhe (rechts/links)	30 x 50 Paar	789.03
M6379.25xx	Sterile Untersuchungshandschuhe (rechts/links)	25 Paar	141.22
M6379.30xx	Sterile Untersuchungshandschuhe (rechts/links)	6x25 Paar	688.20

#### Nitril

Prod.No.	Product	Pack-size	Price in €
M6147.10xx	Alfatex Blau	10 x 100 Stück	78.00
M6146.10xx	Alfatex Natur	10 x 100 Stück	78.00
M6140.10xx	Blossom Aloe Vera Nitril	10 x 100 Stück	106.00
M6148.10xx	Alfatex Blau 30	10 x 100 Stück	88.00

#### Vinyl

Prod.No.	Product	Pack-size	Price in €
M6139.10xx	Blossom Vinyl Handschuh	10 x 100 Stück	82.95
M6139.50xx	Blossom Vinyl Handschuh	50 x 100 Stück	385.88
M6154.10xx	Pretex Vinyl Natur	10 x 100 Stück	51.90
M6141.10xx	Pretex Vinyl Transparent	10 x 100 Stück	48.30

#### Miscellaneous

Prod.No.	Product	Pack-size	Price in €
M6381.12xx	Sterile CoPolymer Einmalhandschuhe	12 x 100 Paar	105.58
M6380.24xx	Sterile CoPolymer Einmalhandschuhe	24 x 100 Stück	149.23
M6377.06xx	Sterile Neopren Operationshandschuhe (rechts/links)	6 x 50 Paar	525.80
M6377.30xx	Sterile Neopren Operationshandschuhe (rechts/links)	30 x 50 Paar	2522.26
M6150.100L	Gummihandschuh	10 Paar	4.84
M6150.100M	Gummihandschuh	10 Paar	4.84
M6150.100S	Gummihandschuh	10 Paar	4.84
M6150.10XL	Gummihandschuh	10 Paar	4.84
M6149.100L	Baumwoll Strickhandschuh.	10 Paar	13.10
M6149.100M	Baumwoll Strickhandschuh.	10 Paar	10.85

## IX. Consumables

M6149.100S	Baumwoll Strickhandschuh.	10 Paar	10.85
M6149.10XL	Baumwoll Strickhandschuh.	10 Paar	22.40

### Protective Masks

Prod.No.	Product	Pack-size	Price in €
M6167.0050	Dust Mask	100 Stück	7.61
M6167.1000	Dust Mask	10 x 100 Stück	70.88
M6166.0100	Face Mask	100 Stück	9.82
M6166.1000	Face Mask	10 x 100 Stück	86.10
M6165.0201	Protective Mask Kim (FPP1)	20 Stück	37.54
M6165.4001	Protective Mask Kim (FPP1)	400 Stück	531.04
M6165.0202	Protective Mask Kim (FPP2)	20 Stück	40.27
M6165.4002	Protective Mask Kim (FPP2)	400 Stück	573.14
M6165.0203	Protective Mask Kim (FPP3)	20 Stück	44.46
M6165.4003	Protective Mask Kim (FPP3)	400 Stück	626.05
M6161.0101	Protective Mask Rob (FFP1)	10 Stück	34.65
M6161.2001	Protective Mask Rob (FFP1)	200 Stück	488.25
M6161.0102	Protective Mask Rob (FFP2)	10 Stück	37.54
M6161.2002	Protective Mask Rob (FFP2)	200 Stück	531.04
M6161.0103	Protective Mask Rob (FFP3)	10 Stück	40.27
M6161.2003	Protective Mask Rob (FFP3)	200 Stück	573.14

### Tubes

#### Cryo Tubes

Prod.No.	Product	Pack-size	Price in €
I2103.000X	Colour code Inserts	5 x 100 pieces	17.96
I2102.1001	Cryotubes 1.2 mL (non-sterile)	1 x 1000 tubes	210.57
I2102.1002	Cryotubes 1.2 mL (non-sterile)	1 x 1000 tubes	216.77
I2102.1000	Cryotubes 1.2 mL (sterile)	10 x 100 tubes	216.77
I2102.2001	Cryotubes 2.0 mL (non-sterile)	1 x 1000 tubes	216.77
I2102.2002	Cryotubes 2.0 mL (non-sterile)	1 x 1000 tubes	216.77
I2102.2000	Cryotubes 2.0 mL (sterile)	10 x 100 tubes	222.96
I2102.4001	Cryotubes 4.0 mL (non-sterile)	1 x 1000 tubes	253.93
I2102.4002	Cryotubes 4.0 mL (non-sterile)	1 x 1000 tubes	253.93
I2102.4000	Cryotubes 4.0 mL (sterile)	5 x 100 tubes	130.06
I2102.5001	Cryotubes 5.0 mL (non-sterile)	1 x 1000 tubes	253.93
I2102.5002	Cryotubes 5.0 mL (non-sterile)	1 x 1000 tubes	253.93
I2102.5000	Cryotubes 5.0 mL (sterile)	5 x 100 tubes	136.25

#### Falcon Type Tubes

Prod.No.	Product	Pack-size	Price in €
I2202.0750	Conical centrifugation tubes (Falcon(TM) type) of 15mL	750 tubes	145.81
I2200.0600	Conical centrifugation tubes (Falcon(TM) type) of 50mL	600 tubes	115.76
I2201.0500	Skirted centrifugation tubes (Falcon(TM) type) of 50mL	500 tubes	105.29
I2203.0750	Sterile conical centrifugation tubes (Falcon(TM) type) of 15mL.	750 tubes	167.37
I2212.0360	Sterile conical centrifugation tubes (Falcon(TM) type) of 50mL	360 tubes	86.63

I2204.0500	Sterile, skirted centrifugation tubes (Falcon(TM) type) of 50mL	500 tubes	123.11
------------	---	-----------	--------

**Reaction Tubes**

Prod.No.	Product	Pack-size	Price in €
I2220.1015	Reaction tube (1.5mL). White	1000 tubes	21.23
I2220.1020	Reaction tube (2.0mL).	1000 tubes	25.95
I2020.0500	Screw Cap Micro Tubes (0.5mL)	500 tubes	28.88
I2021.0500	Screw Cap Micro Tubes (1.5mL)	500 tubes	28.88
I2022.0500	Screw Cap Micro Tubes (2.0mL)	500 tubes	28.88
I2025.0500	Screw caps with O-ring (Blue)	500 caps	52.50
I2024.0500	Screw caps with O-ring (Green)	500 caps	52.50
I2023.0500	Screw caps with O-ring (mixed colours)	500 caps	52.50
I2027.0500	Screw caps with O-ring (Orange)	500 caps	52.50
I2026.0500	Screw caps with O-ring (Red)	500 caps	52.50
I2028.0500	Screw caps with O-ring (Transparent)	500 caps	52.50
I2029.0500	Screw caps with O-ring (Violet)	500 caps	52.50
I2033.0500	Screw caps with O-ring (Yellow)	500 caps	52.50

**PCR Tubes**

Prod.No.	Product	Pack-size	Price in €
I2221.1012	PCR reaction tubes (single tube, domed cap), 0.2 mL	1000 tubes	53.09
I2221.1002	PCR reaction tubes (single tube, flat cap), 0.2 mL	1000 tubes	53.09
I2221.1005	PCR reaction tubes (single tube, flat cap), 0.5 mL	1000 tubes	53.09

**Pipet Tips & Serological Pipets****Pipet Tips**

Prod.No.	Product	Pack-size	Price in €
I2224.2000	0.1 - 10µL pipet tips (non-sterile)	bag of 2000	31.12
I2224.0960	0.1 - 10µL pipet tips (non-sterile)	10 racks	29.77
I2224.9600	0.1 - 10µL pipet tips (non-sterile)	100 racks	275.26
I2225.0960	0.1 - 10µL pipet tips (sterile)	10 racks	33.42
I2225.9600	0.1 - 10µL pipet tips (sterile)	100 racks	308.96
I2226.2000	0.5 - 10µL pipet tips (non-sterile)	bag of 2000	31.12
I2226.0960	0.5 - 10µL pipet tips (non-sterile)	10 racks	29.77
I2226.9600	0.5 - 10µL pipet tips (non-sterile)	100 racks	275.26
I2227.0960	0.5 - 10µL pipet tips (sterile)	10 racks	33.42
I2227.9600	0.5 - 10µL pipet tips (sterile)	100 racks	308.96
I2228.2000	1 - 200µL pipet tips (non-sterile). Natural coloured	bag of 2000	21.07
I2228.0960	1 - 200µL pipet tips (non-sterile). Natural coloured	10 racks	24.16
I2228.9600	1 - 200µL pipet tips (non-sterile). Natural coloured	100 racks	219.08
I2230.2000	1 - 200µL pipet tips (non-sterile). Yellow coloured.	bag of 2000	21.07
I2230.0960	1 - 200µL pipet tips (non-sterile). Yellow coloured.	10 racks	24.16
I2230.9600	1 - 200µL pipet tips (non-sterile). Yellow coloured.	100 racks	219.08
I2229.0960	1 - 200µL pipet tips (sterile). Natural coloured.	10 racks	28.09
I2229.9600	1 - 200µL pipet tips (sterile). Natural coloured.	100 racks	275.26
I2234.1000	100 - 1000µL pipet tips (non-sterile). Blue coloured.	bag of 1000	15.95
I2234.0960	100 - 1000µL pipet tips (non-sterile). Blue coloured.	10 racks	27.81

## IX. Consumables

### Pipet Tips

Prod.No.	Product	Pack-size	Price in €
I2234.9600	100 - 1000µL pipet tips (non-sterile). Blue coloured.	100 racks	235.94
I2232.1000	100 - 1000µL pipet tips (non-sterile). Natural coloured.	bag of 1000	17.70
I2232.0960	100 - 1000µL pipet tips (non-sterile). Natural coloured.	10 racks	30.90
I2232.9600	100 - 1000µL pipet tips (non-sterile). Natural coloured.	100 racks	287.34
I2233.0960	100 - 1000µL pipet tips (sterile). Natural coloured.	10 racks	36.23
I2233.9600	100 - 1000µL pipet tips (sterile). Natural coloured.	100 racks	331.99

### Filter Tips

Prod.No.	Product	Pack-size	Price in €
I2236.5000	Filter Tips (0.1 - 10µL). Sterile.	5000 tips	190.47
I2237.5000	Filter Tips (0.1 - 20µL). Sterile.	5000 tips	190.47
I2239.5000	Filter Tips (1 - 100µL). Sterile.	5000 tips	200.21
I2240.5000	Filter Tips (1 - 200µL). Sterile.	5000 tips	200.21
I2238.5000	Filter Tips (1 - 50µL). Sterile.	5000 tips	200.21
I2241.5000	Filter Tips (100 - 1000µL). Sterile.	5000 tips	219.64
I2236.0960	Filter Tips in rack (0.1 - 10µL). Sterile.	10 racks	53.37
I2237.0960	Filter Tips in rack (0.1 - 20µL). Sterile.	10 racks	53.37
I2239.0960	Filter Tips in rack (1 - 100µL). Sterile.	10 racks	53.37
I2240.0960	Filter Tips in rack (1 - 200µL). Sterile.	10 racks	53.37
I2238.0960	Filter Tips in rack (1 - 50µL). Sterile.	10 racks	53.37
I2241.0960	Filter Tips in rack (100 - 1000µL). Sterile.	10 racks	57.86

### Dispenser Tips

Prod.No.	Product	Pack-size	Price in €
I2101.0100	0.5 mL Dispenser tips	100 pieces	53.26
I2102.0100	1.25 mL Dispenser tips	100 pieces	53.26
I2105.0100	12.5 mL Dispenser tips	100 pieces	53.26
I2103.0100	2.5 mL Dispenser tips	100 pieces	53.26
I2106.0050	25 mL Dispenser tips	50 pieces	30.35
I2104.0100	5.0 mL Dispenser tips	100 pieces	53.26
I2107.0025	50 mL Dispenser tips	25 pieces	16.10
I2108.0000	sterile Dispenser tips	please request	on request

### Serological Pipets

Prod.No.	Product	Pack-size	Price in €
I2215.1000	Disposable serological pipet (1 mL).	box of 1000	100.33
I2218.0800	Disposable serological pipet (10 mL).	box of 800	207.79
I2216.1000	Disposable serological pipet (2 mL).	box of 1000	110.24
I2219.0600	Disposable serological pipet (25 mL).	box of 400	254.11
I2217.1000	Disposable serological pipet (5 mL).	box of 1000	223.27
I2213.3001	Pasteur pipets (3 mL).	500 pieces	81.38

## Racks & Boxes

### Racks

Prod.No.	Product	Pack-size	Price in €
I2104.0001	Cryotube rack (4x10 positions)	1 rack	22.05
I2104.0012	Cryotube rack (4x10 positions)	12 racks	167.22
I2109.0010	Tube-stand LS-10	1 rack	44.83
I2110.0048	Tube-stand LS-48	1 rack	47.19
I2109.0006	Tube-stand LS-6	1 rack	41.29
I2110.0096	Tube-stand LS-96	1 rack	53.09

### Freezer Boxes

Prod.No.	Product	Pack-size	Price in €
I2015.0002	Freezer Box with attached lid (10 x 10). Blue.	1 box	8.66
I2015.0001	Freezer Box with attached lid (10 x 10). Yellow.	1 box	8.66
I2009.0006	Freezer Box with coding (A1 - I9) - black	1 box	4.15
I2009.0002	Freezer Box with coding (A1 - I9) - blue	1 box	3.10
I2009.0003	Freezer Box with coding (A1 - I9) - green	1 box	3.10
I2009.0005	Freezer Box with coding (A1 - I9) - orange	1 box	3.10
I2009.0004	Freezer Box with coding (A1 - I9) - red	1 box	3.10
I2009.0000	Freezer Box with coding (A1 - I9) - white	1 box	3.10
I2009.0001	Freezer Box with coding (A1 - I9) - yellow	1 box	3.10
I2018.0002	Freezer Box without coding - blue	1 box	2.78
I2018.0003	Freezer Box without coding - green	1 box	2.78
I2018.0005	Freezer Box without coding - orange	1 box	2.78
I2018.0004	Freezer Box without coding - red	1 box	2.78
I2018.0000	Freezer Box without coding - white	1 box	2.78
I2018.0001	Freezer Box without coding - yellow	1 box	2.78
I2005.0002	Storage Box (0.2 mL PCR). Neon blue	1 box	7.11
I2005.0003	Storage Box (0.2 mL PCR). Neon green	1 box	7.11
I2005.0005	Storage Box (0.2 mL PCR). Neon orange.	1 box	7.11
I2005.0004	Storage Box (0.2 mL PCR). Neon red	1 box	7.11
I2005.0001	Storage Box (0.2 mL PCR). Neon yellow.	1 box	7.11
I2015.0003	Storage Box with attached lid (10 x 10). Green.	1 box	8.66
I2015.0005	Storage Box with attached lid (10 x 10). Orange.	1 box	8.66
I2015.0004	Storage Box with attached lid (10 x 10). Red.	1 box	8.66

### Card-board Freezer Boxes

Prod.No.	Product	Pack-size	Price in €
I2006.0007	Card-board box with optimised waterproof coating. 133x133 mm. Blue	1 box	4.53
I2006.0008	Card-board box with optimised waterproof coating. 133x133 mm. Green	1 box	4.53
I2006.0009	Card-board box with optimised waterproof coating. 133x133 mm. Red	1 box	4.53
I2006.0010	Card-board box with optimised waterproof coating. 133x133 mm. White	1 box	2.94
I2006.0006	Card-board box with optimised waterproof coating. 133x133 mm. Yellow	1 box	4.53
I2007.0007	Card-board box with optimised waterproof coating. 136x136 mm. Blue	1 box	4.41
I2007.0008	Card-board box with optimised waterproof coating. 136x136 mm. Green	1 box	4.41

## IX. Consumables

### Card-board Freezer Boxes

Prod.No.	Product	Pack-size	Price in €
I2007.0009	Card-board box with optimised waterproof coating. 136x136 mm. Red	1 box	4.41
I2007.0010	Card-board box with optimised waterproof coating. 136x136 mm. White	1 box	2.84
I2007.0006	Card-board box with optimised waterproof coating. 136x136 mm. Yellow	1 box	4.41
I2008.0003	Grid Divider (card-board - 10 x 10 grid - 25 mm)	1 grid	1.89
I2008.0001	Grid Divider (card-board - 10 x 10 grid - 30 mm)	1 grid	1.53
I2008.0004	Grid Divider (card-board - 10 x 10 grid - 40 mm)	1 grid	1.89
I2008.0005	Grid Divider (card-board - 10 x 10 grid - 65 mm)	1 grid	2.24
I2008.0002	Grid Divider (card-board - 9 x 9 grid - 30 mm)	1 grid	1.53
I2006.0002	Storage Box (card-board box). 133x133 mm. Blue	1 box	3.03
I2006.0003	Storage Box (card-board box). 133x133 mm. Green	1 box	3.03
I2006.0004	Storage Box (card-board box). 133x133 mm. Red	1 box	3.03
I2006.0005	Storage Box (card-board box). 133x133 mm. White	1 box	1.89
I2006.0001	Storage Box (card-board box). 133x133 mm. Yellow	1 box	3.03
I2007.0002	Storage Box (card-board box). 136x136 mm. Blue	1 box	2.84
I2007.0003	Storage Box (card-board box). 136x136 mm. Green	1 box	2.84
I2007.0004	Storage Box (card-board box). 136x136 mm. Red	1 box	2.84
I2007.0005	Storage Box (card-board box). 136x136 mm. White	1 box	1.89
I2007.0001	Storage Box (card-board box). 136x136 mm. Yellow	1 box	2.84

### Tissue Culture

Prod.No.	Product	Pack-size	Price in €
I2210.0100	Cell scrapers (24 cm)	100 scrapers	138.92
I2211.0100	Cell scrapers (30 cm)	100 scrapers	161.03
I2242.1620	PS Petri Dishes (d = 55 mm)	81 x 20 dishes	121.39
I2243.0480	PS Petri Dishes (d = 92 mm)	24 x 20 dishes	47.69
I2246.0240	Sterile cell culture petri dishes (d = 100 mm)	240 dishes	98.47
I2244.0900	Sterile cell culture petri dishes (d = 35 mm)	900 dishes	180.22
I2245.0840	Sterile cell culture petri dishes (d = 60 mm)	840 dishes	185.18

## Pipet Tips & Serological Pipets

### High Quality Filter Tips



The very latest micropore technology is used to ensure that all aerosols and liquids are prevented from entering the pipet. This technology eliminates the need for additives, found in self-sealing type filters (which also have the disadvantage of trapping samples in the tip once the filter has sealed). Genaxxon bioscience filter tips use virgin PE filters without additives (which would be a potential source of contamination), and their micropore technology ensures that the filters will never trap samples in the tip. Should the filter become wet due to inadvertent over-pipetting, it will act as a barrier to the passage of liquid, yet the sample in the tip may still be dispensed in the normal way, and can be re-used without fear of contamination.

**Applications** --- The Genaxxon bioscience Filter Tips are engineered for a wide variety of pipettors and are ideal for laboratories with multiple brands of pipettors - even multi-channel pipettors.

### Specifications

Made from PP. Filter made from PE.

All Genaxxon bioscience filter tips are metal free and ideal for DNA amplification (PCR) and microbiological applications including the handling of radioactive, bio-hazardous and corrosive samples.

## Racks & Boxes

### Cryo Boxes with Grid



Cryo Box made from ultra-low temperature polypropylene for the storage of up to 81 cryovials or 1.5ml/2.0ml microcentrifuge tubes. No acetate grid required as box has moulded grid reference (1-9 and A-I) for ease of sample location. Fixed dividers (approx. 13 mm square per division). Translucent colour lid and frosted writing panels. Ultrafreezable for storage down to -196°C, Fully autoclavable at +121°C. Available in 5 neon colours, natural or black. External Dimensions: 130 x 130 mm square and 50 mm fixed height...

### Specifications

Freezable down to - 196°C

Autoclavable at + 121 °C

Dimensions: 130mm x 130 mm x 50 mm (B x L x H)